

CERTIFICATE OF ANALYSIS
Product # HUMANABSRMC-HI-X1
Human AB Serum Male Pooled Not Filtered – LOT 25010097

<u>CHEMISTRY</u>	<u>RESULTS</u>
Glucose	95 mg/dl
Sodium	128 meq/l
Potassium	3.9 meq/l
Chloride	93 meq/l
BUN	13 mg/dl
Creatinine	0.90 mg/dl
Uric Acid	4.3 mg/dl
Calcium	7.3 mg/dl
Phosphorus	4.2 mg/dl
Total Protein	5.7 g/dl
Albumin	3.3 g/dl
Total Bilirubin	<0.2 mg/dl
ALT (SGPT)	5 IU/L
Alkaline Phosphatase	63 IU/L
AST (SGOT)	<5 IU/L
Cholesterol	140 mg/dl
Triglycerides	124 mg/dl

Expiration Date: 31Jan2030

Each donor unit used in the preparation of this material was tested and found negative for HBsAG, HBV-NAT, HIV-1/2, anti-HCV, HCV-NAT, HIV-1 NAT, and a test for Syphilis per FDA regulations.

Male donors blood type AB, collected in FDA registered commercial donor establishments located in the USA and processed in an FDA registered facility located in the USA, in accordance with applicable regulatory and ethical standards.

This product is being sold for further manufacturing only. It is not to be used in humans or animals. For in vitro use only. The user assumes all responsibility for its usage and disposal in accordance with all regulations. The material should be handled as if capable of transmitting infectious agents. Please use universal precautions. No test method can provide total assurance against Hepatitis Virus, Human Immunodeficiency Virus, or other infectious agents that are absent. All blood products should be handled at the Bio-Safety Level 2.

Jazzmyn Bower
Quality Control Specialist



DATE: 10Feb2025

BioIVT is committed to sourcing human and animal biological material in accordance with the highest ethical and regulatory standards.