

Veterinary Rapid Tests

FLUORESCENCE | COLLOIDAL GOLD

*FIA for Laboratory-Grade Accuracy, On-Site
CG for Simple On-Site Rapid Testing*



FIA (Fluorescence Immunoassay) Rapids with:

- High sensitivity, fluorescence-based detection
- Digital readout with FIA analyzer for precision
- Ideal for clinics & labs needing reliable data
- Multiple veterinary tests on one instrument

CGI (Colloidal Gold based Immunoassay) Rapids with:

- Easy-to-use, no instrument required
- Rapid results within minutes
- Field-ready for farms and clinics
- Cost-effective screening solutions

Veterinary Rapid Test Portfolio



Fluorescence Immunoassay (FIA) Portfolio

Higher sensitivity • Digital results • FIA
Instrument Compatible

Livestock

- Bovine Pregnancy Rapid Test (Fluorescence)
- Brucella Antibody Rapid Test (Cow, Sheep, Goat – Serum / WB)

Swine

- Multi-3 Porcine Diarrhea Antigen Rapid Test Strip
- Multi-5 Bovine Diarrhea Antigen Rapid Test Kit

Poultry

- Avian Influenza Virus Generic Antigen Rapid Test

Canine

- CRP Rapid Test
- Pancreas-specific Lipase Rapid Test
- Parvovirus Antigen Rapid Test
- Distemper Virus Antigen Rapid Test
- Hepatitis Virus Antigen Rapid Test
- Corona Virus Antigen Rapid Test

Feline

- Feline Triple (Herpesvirus, Calicivirus, Panleukopenia) Antibody Rapid Test
- Panleukopenia Virus Antibody Rapid Test
- Panleukopenia Virus Antigen Rapid Test
- Herpesvirus Antibody Rapid Test
- Calicivirus Antibody Rapid Test

Colloidal Gold Rapid Test Portfolio

Fast visual results • Cost-effective • No instrument
required

Livestock

- Bovine M. tuberculosis Antibody Rapid Test
- Bovine Viral Diarrhea Virus Antibody Rapid Test
- Bovine Viral Diarrhea Virus Antigen Rapid Test
- Brucella Antibody Rapid Test (Cow, Sheep, Goat – WB)
- Bovine Pregnancy Rapid Test (Colloidal gold)
- Bovine Pregnancy Rapid Test (Latex Microspheres)
- Sheep/Goat Pregnancy Rapid Test Kit (Latex Microspheres)

Swine

- Multi-3 Porcine Diarrhea Antigen Test Strip

Canine

- Parvovirus Antigen Rapid Test
- Distemper Virus Antigen Rapid Test
- Corona Virus Antigen Rapid Test
- Distemper Virus Antibody Rapid Test
- Hepatitis Virus Antibody Rapid Test

Feline

- Panleukopenia Virus Antibody Rapid Test
- Calicivirus Antibody Rapid Test
- Herpesvirus Antibody Rapid Test
- Panleukopenia Virus Antigen Rapid Test
- Calicivirus Antigen Rapid Test
- Influenza Virus Antigen Rapid Test

From livestock to pets, choose between digital FIA precision or affordable colloidal gold rapid tests — all from one trusted manufacturer.

Herd health, clinical screening, and early diagnosis made easy.