

INVITROGRO™ HI MEDIUM

Product No.	Description	Size
Z99009	INVITROGRO HI Medium, serum-free	250 mL
Z990012	INVITROGRO HI Medium, serum-free	500 mL
Z990002	INVITROGRO HI Medium, serum-free	1 L

Product Description:

INVITROGRO HI Medium is a hepatocyte incubation medium used for culturing primary hepatocytes¹. INVITROGRO HI Medium is serum-free and formulated without antibiotics to increase its shelf life. BioIVT recommends the use of the TORPEDO™ Antibiotic Mix with INVITROGRO HI Medium.

Stability: 6 months for INVITROGRO HI Medium, serum-free, without TORPEDO Antibiotic Mix
7 days for Complete INVITROGRO HI Medium

Storage: 4°C

Related Products:

Product No.	Size
Z99000	TORPEDO Antibiotic Mix 5.5 ml
Z990007	TORPEDO Antibiotic Mix 11 mL
Z990008	TORPEDO Antibiotic Mix 22 mL

Reference:

1. Roymans, D.; Van Looveren, C.; Leone, A.; Parker, J. B.; McMillan, M.; Johnson, M. D.; Koganti, A.; Gilissen, R.; Silber, P.; Mannens, G.; Meuldermans, W. Determination of cytochrome P450 1A2 and cytochrome P450 3A4 induction in cryopreserved human hepatocytes. *Biochem. Pharmacol.* **2004**, 67(3), 427–437.

Caution: This product is being sold for research and/or manufacturing purposes only. The biological samples supplied by BioIVT, or any material isolated from the samples, are for in-vitro research use only and are not to be used as a source of material for clinical therapies. Human material may be used in vivo in animals. The user assumes all responsibility for its usage and disposal, in accordance with all regulations.