

INVITROGRO™ CP NPC MEDIUM

Product No.	Description	Size
Z990031	INVITROGRO CP NPC Medium	250 mL

Product Description:

INVITROGRO CP NPC Medium is used for thawing and plating primary human Kupffer cells. INVITROGRO CP NPC Medium contains serum and is formulated without antibiotics to increase its shelf life. BioIVT recommends the use of the TORPEDO™ Antibiotic Mix with INVITROGRO CP NPC Medium.

Stability: 6 months for INVITROGRO CP NPC Medium without TORPEDO Antibiotic Mix
7 days for Complete INVITROGRO CP NPC Medium

Storage: 2-8°C

Related Products:

Product No.	Size
Z99000	TORPEDO Antibiotic Mix 5.5 ml

Caution: This product is being sold for research and/or manufacturing purposes only. The biological samples supplied by BioIVT, or any material isolated from the samples, are for in-vitro research use only and are not to be used as a source of material for clinical therapies. Human material may be used in vivo in animals. The user assumes all responsibility for its usage and disposal, in accordance with all regulations.