

INVITROGRO™ HT MEDIUM

Product No.	Description	Size
Z99019	INVITROGRO HT Medium	250 mL
Z990005	INVITROGRO HT Medium	500 mL
Z990006	INVITROGRO HT Medium	1 L

Product Description:

INVITROGRO HT Medium is a hepatocyte incubation medium used for culturing primary hepatocytes¹. INVITROGRO HT Medium contains serum and is formulated without antibiotics to increase its shelf life. BioIVT recommends the use of the TORPEDO™ Antibiotic Mix with INVITROGRO HT Medium.

Stability: 6 months for INVITROGRO HT Medium without TORPEDO Antibiotic Mix
7 days for Complete INVITROGRO HT Medium

Storage: 4°C

Related Products:

Product No.	Size
Z99000	TORPEDO Antibiotic Mix 5.5 ml
Z990007	TORPEDO Antibiotic Mix 11 mL
Z990008	TORPEDO Antibiotic Mix 22 mL

Reference:

1. Prabhu, S.; Fackett, A.; Lloyd, S.; McClellan, H. A.; Terrell, C. M.; Silber, P. M.; Li, A. P. Identification of glutathione conjugates of troglitazone in human hepatocytes. *Chem. Biol. Interact.* **2002**, 142(1–2), 83–97.

Caution: This product is being sold for research and/or manufacturing purposes only. The biological samples supplied by BioIVT, or any material isolated from the samples, are for in-vitro research use only and are not to be used as a source of material for clinical therapies. Human material may be used in vivo in animals. The user assumes all responsibility for its usage and disposal, in accordance with all regulations.